

**TANOPLUS**  
AR-fresh *by* partoeno

*The Art  
and Science  
of oenology*

**TANOPLUS**  
AR-fresh *by* partoeno

## FORMULATION OF PROANTHOCYANIDIC PLANT EXTRACT AND EUCALYPTUS GLOBULUS FOR HUMAN CONSUMPTION

**STABILIZES THE COLOR. REDEVELOPING THE FRESHNESS OF HEAVY WINES.**

100% NATURAL, WITHOUT PRESERVATIVES. WITHOUT GMO.



### Application

TANOPLUS AR FRESH® contributes to the stabilization of the color.

TANOPLUS AR FRESH® is advised for heavy beverages and/or destructured wines.

### Description and Uses

TANOPLUS AR FRESH® is a product developed by Partoeno R&D and tested by the University of Bordeaux, France and key wine consultants on different wines, regions and vintages.

#### • Another new product ?

Partoeno has compared TANOPLUS AR FRESH® to most of existing formulations on the market. In all cases, it was judged by an independent panel as giving better results and being more cost-effective.

**TANOPLUS AR FRESH® IS A PRODUCT INTENDED TO REFRESH THE WINES AND ADVICE FOR HEAVY WINES.**

### When to use TANOPLUS AR FRESH® ?

Wines that are considered to be outside the acceptable range eg : all wines lacking of structure should be lab- treated with TANOPLUS AR FRESH®.

The treated wines can be tasted the same day of the treatment, however it is recommended to wait at least 48 hours.

Furthermore, revalue wines of average range at a lower cost.

It is recommended to use TANOPLUS AR FRESH® at a dosage between **2 et 15 g/hl** for red wine and between **2 et 7 g/hl** for white and rosé wine after dilution in wine. The dosage, and therefore the treatment cost, can be optimised by doing direct trials in bottles at different concentrations.

**TANOPLUS**  
AR-fresh *by* partoeno

### Characteristic physico-chemical and limit bacteriological

Aspect :	fine powder	Quality of polyphenols :	min. 65%
Color :	Orange coloured	Mercury:	< 1 ppm
Smell :	characteristic of Eucalyptus	Leads :	< 5 ppm
Arsenic :	< 3 ppm		
Humidity :	< 8 %		
Toxines et mycotoxines :	abs.	Escherichia coli :	abs.
Salmonella :	abs.	Coliforms :	abs.
Staphylococcus aureus :	abs.		

### Shelf life, Storage and Packaging

3 years ( from date of production ) in their original packaging, stored in a cool, odourless and dry place.

Pack de 25kgs. ( Other packaging, ask us ).

For more information, ask the Material Safety Data Sheet TANOPLUS AR FRESH®  
CTN : 32011000