

FLAVTAN
by
VQR *partoeno*

*The Art
and Science
of oenology*

FLAVTAN
by
VQR *partoeno*

PURE PROANTHOCYANIDIC TANNIN :

Quebracho extract for human consumption.

FLAVTAN VQR® CONTRIBUTES TO THE STABILIZATION OF COLOUR.

100% NATURAL, WITHOUT PRESERVATIVES. WITHOUT GMO.



Application

This tannin is issued from the Quebracho Sinopsis Lorentzii (South America). FLAVTAN VQR® contributes to the clarification and the stabilization of wines.

Description and Uses

FLAVTAN VQR® is a product developed by Partoeno R&D and tested by the University of Bordeaux, France and key wine consultants on different wines, regions and vintages. This product is very reactive to the proteins.

This product stabilizes the potential of oxidoreduction.

FLAVTAN VQR® IS USED FOR WINES PRESENTING DIFFICULTIES OF CLARIFICATION AND CONTRIBUTES TO THE STABILITY OF COLOUR.

When to use FLAVTAN VQR®?

Wines that are considered to be outside the acceptable range eg : all wines presenting difficulties of

clarification and stabilization. The treated wines can be tasted the same day of the treatment, however it is recommended to wait at least 48 hours.

Furthermore, revalue wines of average range at a lower cost.

FLAVTAN VQR® have a very important anti-oxidation power.

It is recommended to use FLAVTAN VQR® at a dosage between **2 to 30 g/hl** of red and rosé wine after dilution in wine. The dosage, and therefore the treatment cost, can be optimised by doing direct trials in bottles at different concentrations.

FLAVTAN
by
VQR *partoeno*

Characteristic physico-chemical and limit bacteriological

Aspect :	fine powder	Quality of polyphenols :	min. 65%
Color :	orange/red coloured	Arsenic :	< 3 ppm
Mercury :	< 1ppm	Leads :	< 5 ppm
Humidity :	< 8 %	Iron :	< 50 ppm
Toxines et mycotoxines :	abs.	Escherichia coli :	abs.
Salmonella :	abs.	Coliforms :	abs.
Staphylococcus aureus :	abs.	Extract Solvent :	Alcohol

Shelf life, Storage and Packaging

3 years (from date of production) in their original packaging, stored in a cool, odourless and dry place.

Pack de 25kgs. (Other packaging, ask us).

For more information, ask the Material Safety Data Sheet FLAVTAN VQR®
CTN : 32011000